

AUTOMATED MICROSCOPIC IMAGE ACQUISITION, COMPOSITING, AND DISPLAY

ABSTRACT

5 An automated microscope and computer system captures a set of images for
a capture area in a plurality of focal planes. The images can then be integrated into
composite images for browsing to simulate viewing an item, such as a biological
sample, under a microscope. A corrective filter can be constructed from the images
to avoid an effect called "tiling." Before capture, variable focal plane error can be
10 avoided by collecting z locations for a set of points in the capture area. During
image browsing, entire composite images can be loaded into memory in compressed
form. Compressed image portions can be pre-decompressed to avoid delay as a
browsing user navigates throughout the composite images. Pre-decompression can
be done by a thread separate from the thread performing navigation operations.

15